


FREE-ANGLE BIOPSY GUIDES

**Thyroid and breast interventional
disposable ultrasound biopsy guides**

Used for thyroid FNA, thyroid biopsy, and breast biopsy.

- The free-angle ensures complete freedom during the procedure, ensuring a precise display of the needle on the ultrasound image, making superficial interventions more accurate and safe.
- 6 kinds of needle guides with different specifications, for various clinical applications such as thyroid FNA, thyroid and breast biopsy, and treatment.
- Designed for a one-time use with a self-destruction function, in order to reduce the risk of infection and cross-contamination.



Contact us

E-mail: sales@leapmed.com
Fax: +86 0512 6655 4498-8002
Address: Wuzhong Science Park, building 2
North Guandu Road no.38
Wuzhong District, Suzhou, China
Phone: **+86 0512 66554498**



Official website



Facebook page

Step by step user guide

01.



Install the probe cover

The ultrasound coupling agent is smeared on the probe, which is then covered by the disposable probe cover.


02.



Install the biopsy guide

Place the biopsy guide on the ultrasound probe, and push the lock to stabilize it.

03.



Select the proper needle adapter

Select the adapter that matches the diameter of your needle.

04.



During the surgery, the guide can separate with the needle

At any point during the surgery, the angle of the needle can be changed, and the guide can be taken away from the needle.

05.



Remove the guide

Pull the lock and dispose of the biopsy guide.

Those steps are described for reference only, they do not prevail the instructions of any medical procedure.

More disposable products are currently under development, if you can't find a guide that meets your needs within the disposable biopsy guide series, please note that our metal biopsy guides can be customized.

Compatible probes

GE Healthcare



GE 9L/9L-D GE ML6-15-D



PHILIPS

Philips L12-5
Philips CX L12-3/CV L12-4

SIEMENS

Siemens 9L4/VF10-5

**HITACHI
ALOKA**

Hitachi L74M/L55
Hitachi L75/L64

esaote
Creativity in Healthcare

Esaote LA523 Esaote LA332

mindray
healthcare within reach

Mindray L12-3E/7L4A/7L4

SonoSite
Imagine What You Can Do

Sonosite HFL38X
Sonosite C60/C60e/C60X